



## Specifications: So-Lo & So-Lo with Shower Deck Models

<b>Capacity</b>	1 Adult
<b>Dimensions</b>	<b>So-Lo:</b> 32"W x 59 5/8"L x 18"H <b>So-Lo with Shower Deck:</b> 44"W x 59.5"L x 18"H
<b>Door Swing</b>	Right or Left Hand
<b>Warranty</b>	<p>*Lifetime warranty on door seal – first year parts &amp; labor, remaining years, parts only.</p> <p>*3 year warranty on bathtub, door and latching mechanism- 1 year parts and labor, remaining 2 years parts only.</p> <p>*One year warranty on hydrotherapy components, heated seat, flip seat and latch, and drain components purchased from Safety Bath</p> <p>*Warranty does not cover the installation of your So-Lo bathtub, nor does it cover modifications made to the flip-down seat.</p>
<b>Length</b>	59 5/8" L
<b>Width</b>	<b>So-Lo:</b> 32" <b>So-Lo with Shower Deck:</b> 44"
<b>Height</b>	<b>So-Lo:</b> 18" <b>So-Lo with Shower Deck:</b> 18"
<b>Step in height</b>	4"
<b>Door opening width</b>	14.25"
<b>Weight</b>	Without Jets: 90 lbs With jets – 110 lbs
<b>Door seal</b>	Extruded silicone rubber gasket
<b>So-Lo Safety Seat</b>	<p>*Measures: 18"W x 30"L</p> <p>*Comes with a factory mounted anchor system and 4 screw caps</p> <p>* Safety Seat can be mounted after the So-Lo is installed</p> <p>*Flip-down handle locks seat in upright position when not in use</p> <p>*Constructed of HDPE High Density Polyethylene</p>
<b>Tub surface</b>	Acrylic - "Kohler" white

<b>Tub construction</b>	Acrylic with fiberglass reinforcement
<b>Hinges</b>	Stainless steel
<b>Drain size</b>	1.5", ABS connection
<b>Drain Location</b>	Right or Left hand drain available
<b>Access panel</b>	Measures: 8" x 8" Located in front <b>on jetted tub only</b> , opposite drain location
<b>Heated Seat</b>	<b>So-Lo with Shower Deck</b> model only * Shower deck measures 12"W x 59 5/8"L * 110/120 V – 15 A GFI outlet required
<b>Power Input</b>	110/120 V – 15 A GFI outlet required for jetted tubs
<b>Piping</b>	PVC flex line
<b>Jets</b>	<b>Option:</b> * 12 Warm Air Jets - 12 chrome jets with one way check valves & built-in purge cycle after tub is completely drained.